


INSTRUCTIONS FOR RETRIEVING SAVER™ EVO FILE


To retrieve a Saver™ EVO file from the samaritan® PAD:


1. Connect the samaritan® PAD to the PC using the supplied data cable.
2. Open the Saver™ EVO software on your PC. (If you do not have Saver™ EVO software installed on your PC, download the software from <http://uk.heartsine.com/support/upload-saver-evo/>)


3. Click  button to connect to the samaritan PAD.

4. Click  button to display a list of saved events.

5. From the list, highlight the event you want to download.

6. Click  button to view the event.


7. Click  button to save the file and enter patient information:




The screenshot shows the 'Enter Patient Details' dialog box in the Saver™ EVO software. The dialog box is overlaid on a background window showing a waveform and status indicators. The dialog box contains the following fields and buttons:


- Record Type: Subject (dropdown menu)
- Surname: John (text input)
- Forename: Smith (text input)
- D.O.B.: 01/01/44 (text input)
- Patient ID: 12345 (text input)
- Event Date and Time: 01 Nov 2012 14:20:15 GMT Standard Time (text input)
- Buttons: New Record Type, Save, Cancel

8. After you enter the patient detail, click  button to finish saving the event to the PC.

9. Click  button three times to return to the Home screen.

10. Click  button to open the Saved Events directory.

11. Select the required event file.

12. Click  button to export the event file.

13. Enter a file name and save the Saver™ EVO file to an accessible location on the PC.

14. Email the Saver™ EVO file, along with the Customer Event Form, to uploads@heartsine.com.

Complete Saver™ EVO manuals are available at <http://www.heartsine.com/support/product-manuals/>.

